

## Low level fume cupboard (DIGIM 12™ NP, DIGIM 15™ NP, DIGIM 18™ NP)




### → GENERAL PROPERTIES

- Made of galvanized metal sheet epoxy powder coated (higher resistance)
- Completely in accordance with HRN EN 14725-2
- Different types of worktops (technical ceramic material, epoxy, polypropylene, stainless steel),
- Safety non-dispersing 6mm glass
- Supply fittings (water, gas, vacuum- optional, demineralised water-optional, oxygen-optional, etc.)
- 220V16A socket (standard 3 pieces) IP44 protection
- Ventilator with engine in EX version-optional (Exe II T1-T3 protection, IP54, Zone 1 and 2)
- AFM Processor management- optional
- According to HRN D.E2.105. EN 14175-2

### → DIMENSIONS ACCORDING TO EN norms 14725-2

	DIGIM 12™ NP	DIGIM 15™ NP	DIGIM 18™ NP
<b>width</b>	1200	1500	1800
<b>hight</b>	2500	2500	2500
<b>tot.hight<sup>1</sup></b>	2800	2800	2800
<b>hight wt<sup>2</sup></b>	500	500	500
<b>width wt<sup>3</sup></b>	1000	1300	1600
<b>depth</b>	900(950)	900(950)	900(950)

1 - necessary height, 2 -worktop height above the floor, 3 - working width



## Low level fume cupboard (DIGIM 12™ NP, DIGIM 15™ NP, DIGIM 18™ NP)

### → TECHNICAL CERAMIC MATERIAL

- Worktop and tub
- High chemical and biological resistance
- No joints or junctions
- Elevated edges on all four sides

### → SAFETY CUPBOARDS

- For acids, bases and inflammable substances
- Self-ventilated
- Top quality
- TUV certificate

### → PROCESSOR MANAGEMENT

- Complete fume cupboard work control
- Air circulation control according to EN 12924 II part
- Optimization of ventilator rotation - saving heating at workplace

### → CARBON FILTERS

- different dimensions and purposes
- high ventilation efficiency
- simple and efficient maintenance



### Grga & Melita

Croatia, 10 000 Zagreb, Sveti Duh 62

tel: +385 (1) 3745 510, 3906 616

mob: +385 (98) 355 274, 492 717

fax: +385 (1) 3906 615

e-mail: [grga-melita@zg.t-com.hr](mailto:grga-melita@zg.t-com.hr)

web: [www.grga-melita.hr](http://www.grga-melita.hr)