

Laboratory cupboards



GENERAL SPECIFICATIONS

The cupboard body is made of moisture-proof RAL 7047 chipboard.

All edges are protected by 0.8 mm ABS tape. The moisture-proof MDF rear panel is also coated with RAL 7047 melamine foil on both sides.

Doors are made of moisture-proof chipboard or framed. All drawers are made of metal and equipped with soft-closing system. 150 mm cupboard levelling legs are made of RAL 7047 polypropylene.

Ergonomic handles are made from ABS with labelling aperture.

